


Product Datasheet

Porcine F1+2 ELISA Kit (orb566493)

Description	Porcine F1+2 ELISA Kit
Reactivity	Porcine
Range	6.25-400ng/ml
Target	F1+2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	3.75ng/ml
Expiration Date	Please enquire.



Porcine
F1+2(Prothrombin
Fragment 1+2)
E...