


Product Datasheet

Porcine FSH ELISA Kit (orb566503)

Description	Porcine FSH ELISA Kit
Reactivity	Porcine
Range	23.438-1500ng/ml
Target	FSH
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01228
Sensitivity	14.063ng/ml
Expiration Date	Please enquire.



Porcine FSH(Follicle-Stimulating Hormone...