


Product Datasheet

Porcine GC ELISA Kit (orb566504)

Description	Porcine GC ELISA Kit
Reactivity	Porcine
Range	62.5-4000pg/ml
Target	GC
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01274
Sensitivity	37.5pg/ml
Expiration Date	Please enquire.



Porcine
GC(Glucagon)
ELISA Kit