


Product Datasheet

Porcine GH ELISA Kit (orb566507)

Description	Porcine GH ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/ml
Target	GH
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01248
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Porcine
GH(Growth
Hormone)
ELISA Kit