


Product Datasheet

Porcine HIF-1alpha ELISA Kit (orb566510)

Description	Porcine HIF-1alpha(Hypoxia Inducible Factor 1 Alpha) ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/ml
Target	HIF-1 α
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	A0T1H6
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Porcine HIF-1 α (Hypoxia Inducible Factor...)