


Product Datasheet

Porcine IFN-gamma ELISA Kit (orb566515)

Description	Porcine IFN-gamma(Interferon Gamma) ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/ml
Target	IFN- γ
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P17803
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Porcine IFN- γ (Interferon Gamma) ELISA K...