


Product Datasheet

Porcine IL-10 ELISA Kit (orb566526)

Description	Porcine IL-10 ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/ml
Target	IL-10
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q29055
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Porcine IL-10(Interleukin-10) ELISA Kit