


Product Datasheet

Porcine IL-17 ELISA Kit (orb566530)

Description	Porcine IL-17 ELISA Kit
Reactivity	Porcine
Range	15.625-1000pg/ml
Target	IL-17
Storage	2 - 8 °C for 12 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q60I29
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Porcine IL-17(Interleukin 17) ELISA Kit

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)