


## Product Datasheet

### Porcine IL-18 ELISA Kit (orb566531)

<b>Description</b>	Porcine IL-18 ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	IL-18
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>O19073</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Porcine IL-18(Interleukin 18) ELISA Kit