


## Product Datasheet

### Porcine IL-1alpha ELISA Kit (orb566533)

<b>Description</b>	Porcine IL-1alpha(Interleukin 1 Alpha) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	IL-1 $\alpha$
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>P18430</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Porcine IL-1 $\alpha$ (Interleukin 1 Alpha) ELIS...