


Product Datasheet

Porcine IL-1beta ELISA Kit (orb566534)

Description	Porcine IL-1beta ELISA Kit
Reactivity	Porcine
Range	62.5-4000pg/ml
Target	IL-1 β
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P26889
Sensitivity	37.5pg/ml
Expiration Date	Please enquire.



Porcine IL-1 β (Interleukin 1 Beta) ELISA...