


Product Datasheet

Porcine IL-4 ELISA Kit (orb566537)

Description	Porcine IL-4 ELISA Kit
Reactivity	Porcine
Range	15.625-1000pg/ml
Target	IL-4
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q04745
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Porcine IL-4(Interleukin 4) ELISA Kit