


## Product Datasheet

### Porcine LH ELISA Kit (orb566541)

<b>Description</b>	Porcine LH ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	0.781-50ng/ml
<b>Target</b>	LH
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>P01232</b>
<b>Sensitivity</b>	0.469ng/ml
<b>Expiration Date</b>	Please enquire.



Porcine  
LH(Luteinizing  
Hormone)  
ELISA Ki...