


Product Datasheet

Porcine MIP-1beta ELISA Kit (orb566546)

Description	Porcine MIP-1beta ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/ml
Target	MIP-1 β
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q711P4
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Porcine MIP-1 β (Macrophage Inflammatory

...