


## Product Datasheet

### Porcine NOS3 ELISA Kit (orb566554)

<b>Description</b>	Porcine NOS3 ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	12.5-800pg/ml
<b>Target</b>	NOS3/eNOS
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>Q28969</b>
<b>Sensitivity</b>	7.5pg/ml
<b>Expiration Date</b>	Please enquire.



Porcine NOS3/eNOS(Nitric Oxide Synthase ...