


Product Datasheet

Porcine OPN ELISA Kit (orb566561)

Description	Porcine OPN ELISA Kit
Reactivity	Porcine
Range	62.5-4000pg/ml
Target	OPN
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P14287
Sensitivity	37.5pg/ml
Expiration Date	Please enquire.



Porcine
OPN(Osteopontin)
ELISA Kit