


Product Datasheet

Porcine PDGFB ELISA Kit (orb566564)

Description	Porcine PDGFB ELISA Kit
Reactivity	Porcine
Range	46.875-3000pg/ml
Target	PDGFB
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P20034
Sensitivity	28.125pg/ml
Expiration Date	Please enquire.



Porcine PDGFB(Platelet Derived Growth Fa...