


Product Datasheet

Porcine S100A6 ELISA Kit (orb566571)

Description	Porcine S100A6 ELISA Kit
Reactivity	Porcine
Range	62.5-4000pg/ml
Target	S100A6
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q2EN75
Sensitivity	37.5pg/ml
Expiration Date	Please enquire.



Porcine
S100A6(S100
Calcium
Binding
Prot...