


Product Datasheet

Porcine TNF-alpha ELISA Kit (orb566582)

Description	Porcine TNF-alpha(Tumor Necrosis Factor Alpha) ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/ml
Target	TNF- α
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P23563
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Porcine TNF- α (Tumor
Necrosis
Factor Alp...