


Product Datasheet

Porcine VWF ELISA Kit (orb566589)

Description	Porcine VWF ELISA Kit
Reactivity	Porcine
Range	1.563-100ng/ml
Target	VWF
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q28833
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Porcine
VWF(Von
Willebrand
Factor)
ELISA...