


## Product Datasheet

### Porcine CLDN7 ELISA Kit (orb566600)

<b>Description</b>	Porcine CLDN7 ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	125-8000pg/ml
<b>Target</b>	CLDN7
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>H0W9G6</b>
<b>Sensitivity</b>	75pg/ml
<b>Expiration Date</b>	Please enquire.



Porcine  
CLDN7(Claudin-  
7) ELISA Kit