


Product Datasheet

Porcine COR ELISA Kit (orb566639)

Description	Porcine COR ELISA Kit
Reactivity	Porcine
Range	0.391-25ng/ml
Target	COR
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	0.234ng/ml
Expiration Date	Please enquire.



Porcine
COR(Cortisol)
ELISA Kit