


## Product Datasheet

### Porcine HMGB1 ELISA Kit (orb566644)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Porcine HMGB1 ELISA Kit  |
| <b>Reactivity</b>      | Porcine  |
| <b>Range</b>           | 31.25-2000pg/ml  |
| <b>Target</b>          | HMGB1  |
| <b>Storage</b>         | 2-8 °C for 6 months  |
| <b>Note</b>            | For research use only  |
| <b>Sample Types</b>    | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| <b>Uniprot ID</b>      | <b>P12682</b>  |
| <b>Sensitivity</b>     | 18.75pg/ml   |
| <b>Expiration Date</b> | Please enquire.  |



Porcine HMGB1(High mobility group protei...