


Product Datasheet

Rat IP-10 ELISA Kit (orb566674)

Description	Rat IP-10 ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	IP-10/CXCL10
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P48973
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat IP-10/CXCL10(Interferon Gamma Induce...