


Product Datasheet

Rat IGF-1 ELISA Kit (orb566675)

Description	Rat IGF-1 ELISA Kit
Reactivity	Rat
Range	62.5-4000pg/ml
Target	IGF-1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P08025
Sensitivity	37.5pg/ml
Expiration Date	Please enquire.



Rat IGF-1 (Insulin-Like Growth Factor 1) ...