


Product Datasheet

Rat IGF-2 ELISA Kit (orb566676)

Description	Rat IGF-2 ELISA Kit
Reactivity	Rat
Range	46.875-3000pg/ml
Target	IGF-2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01346
Sensitivity	28.125pg/ml
Expiration Date	Please enquire.



Rat IGF-2(Insulin Like Growth Factor 2) ...