


Product Datasheet

Rat IGFBP-5 ELISA Kit (orb566677)

Description	Rat IGFBP-5 ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/ml
Target	IGFBP-5
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P24594
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Rat IGFBP-5 (Insulin Like Growth Factor B...