


Product Datasheet

Rat IL6 ELISA Kit (orb566685)

| | |
|------------------------|--|
| Description | Rat IL6 ELISA Kit |
| Reactivity | Rat |
| Range | 62.5-4000pg/ml |
| Target | IL-6 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Sample Types | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| Uniprot ID | P20607 |
| Sensitivity | 37.5pg/ml |
| Expiration Date | Please enquire. |



Rat IL-6(Interleukin-6) ELISA Kit