


Product Datasheet

Rat MIP-1alpha ELISA Kit (orb566691)

Description	Rat MIP-1alpha/CCL3(Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	MIP-1 α /CCL3
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P50229
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat MIP-1 α /CCL3(Macrophage Inflammatory...)