


Product Datasheet

Rat MIP-1alpha ELISA Kit (orb566691)

Description	Rat MIP-1alpha/CCL3(Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	MIP-1 α /CCL3
Storage	2 - 8 °C for 12 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Assay Time	4 hours
Uniprot ID	P50229
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat MIP-1 α /CCL3(Macrophage Inflammatory Protein 1 Alpha) ELISA Kit

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)