


Product Datasheet

Rat NT-3 ELISA Kit (orb566696)

Description	Rat NT-3 ELISA Kit
Reactivity	Rat
Range	31.2-2000pg/ml
Target	NT-3
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P18280
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat NT-3 (Neurotrophin 3) ELISA Kit