


Product Datasheet

Rat RANTES ELISA Kit (orb566698)

Description	Rat RANTES ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	Rantes
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P50231
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
Rantes(Regulated
On Activation,
Norm...