


Product Datasheet

Rat Tshb ELISA Kit (orb566710)

Description	Rat Tshb ELISA Kit
Reactivity	Rat
Range	0.156-10uIU/ml
Target	Tshb
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P04652
Sensitivity	0.094uIU/ml
Expiration Date	Please enquire.



Rat
Tshb(Thyrotropin
subunit beta)
ELISA...