


Product Datasheet

Rat IL-5 ELISA Kit (orb566723)

Description	Rat IL-5 ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/ml
Target	IL-5
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q08125
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat IL-5 (Interleukin-5) ELISA Kit