


Product Datasheet

Rat PDGF B ELISA Kit (orb566735)

Description	Rat PDGF B ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	PDGFB
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q05028
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat
PDGFB(Platelet-
derived growth
factor...