


Product Datasheet

Rat FGF23 ELISA Kit (orb566737)

Description	Rat FGF23 ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	FGF23
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q8VI82
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
FGF23(Fibroblast
growth factor
23) E...