


Product Datasheet

Rat Neuropeptide Y ELISA Kit (orb566741)

Description	Rat Npy(neuropeptide Y) ELISA Kit
Reactivity	Rat
Range	31.2-2000pg/ml
Target	Npy
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P07808
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat Npy(Pro-neuropeptide Y) ELISA Kit