


Product Datasheet

Rat Gip ELISA Kit (orb566742)

Description	Rat Gip ELISA Kit
Reactivity	Rat
Range	0.625-40ng/ml
Target	Gip
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q06145
Sensitivity	0.375ng/ml
Expiration Date	Please enquire.



Rat
Gip(Gastric
inhibitory
polypeptide)
...