


Product Datasheet

Rat IL24 ELISA Kit (orb566748)

Description	Rat IL24 ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/ml
Target	IL-24
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9JI24
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
IL24(Interleukin-
24) ELISA Kit