


Product Datasheet

Rat PGC ELISA Kit (orb566787)

Description	Rat PGC ELISA Kit
Reactivity	Rat
Range	4.688-300ng/ml
Target	PGC
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P04073
Sensitivity	2.813ng/ml
Expiration Date	Please enquire.



Rat
PGC(Pepsinogen
C) ELISA Kit