


Product Datasheet

Rat AGE ELISA Kit (orb566809)

Description	Rat AGE ELISA Kit
Reactivity	Rat
Range	0.313-20ng/ml
Target	AGE
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P51607
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
AGE(Advanced
glycation end-
product) ...