


Product Datasheet

Rat Reg3b ELISA Kit (orb566834)

Description	Rat Reg3b ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	Reg3b
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P25031
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
Reg3b(Regenerating
islet-derived pro...