


Product Datasheet

Rat Stab2 ELISA Kit (orb566843)

Description	Rat Stab2 ELISA Kit
Reactivity	Rat
Range	0.312-20ng/ml
Target	Stab2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q8CFM6
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
Stab2(Stabilin-
2) ELISA Kit