


Product Datasheet

Rat CYP27B1 ELISA Kit (orb566844)

Description	Rat CYP27B1 ELISA Kit
Reactivity	Rat
Range	0.312-20ng/ml
Target	Vdr
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P13053
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
Vdr(Vitamin
D3 receptor)
ELISA Kit