


Product Datasheet

Rat CD44 ELISA Kit (orb566868)

Description	Rat CD44 ELISA Kit
Reactivity	Rat
Range	0.312-20ng/ml
Target	Cd44
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P26051
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
Cd44(CD44
antigen)
ELISA Kit