


Product Datasheet

Rat Adh1 ELISA Kit (orb566870)

Description	Rat Adh1 ELISA Kit
Reactivity	Rat
Range	0.625-40ng/ml
Target	Adh1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P06757
Sensitivity	0.375ng/ml
Expiration Date	Please enquire.



Rat
Adh1(Alcohol
dehydrogenase
1) ELISA ...