


Product Datasheet

Rat Kik8 ELISA Kit (orb566872)

Description	Rat Kik8 ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	Kik8
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	O88780
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
Kik8(Kallikrein-
8) ELISA Kit