


## Product Datasheet

### Rat Prdx5 ELISA Kit (orb566875)

<b>Description</b>	Rat Prdx5 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.78-50ng/ml
<b>Target</b>	Prdx5
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>Q9R063</b>
<b>Sensitivity</b>	0.469ng/ml
<b>Expiration Date</b>	Please enquire.



Rat  
Prdx5(Peroxiredoxin-  
5, mitochondrial...