


Product Datasheet

Rat Apolipoprotein E ELISA Kit (orb566876)

Description	Rat Apolipoprotein E ELISA Kit
Reactivity	Rat
Range	3.125-200ng/ml
Target	ApoE
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P02650
Sensitivity	1.875ng/ml
Expiration Date	Please enquire.



Rat
ApoE(Apolipoprotein
E) ELISA Kit