


Product Datasheet

Rat TLR4 ELISA Kit (orb566882)

Description	Rat TLR4 ELISA Kit
Reactivity	Rat
Range	0.312-20ng/ml
Target	TLR4
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9QX05
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
TLR4(Toll-
like receptor
4) ELISA Kit